

SURJEET INTERNATIONAL



DEALS IN :
PHYSIOTHERAPY ❖ ELECTRO THERAPY
BEAUTY & SLIMMING ❖ HEALTH & FITNESS

ABOUT US

We **Surjeet International**, are manufacturing, importer, exporter and wholesaling a high-quality Beauty and Physiotherapy Equipment like Laser Hair Removal Machine, Tattoo Removal Machine, Skin Care Machine etc. Offered products range consists of Laser Hair Removal Machine, Tattoo Removal Machine and Skin Care Machine. This range is precisely manufactured using the best quality raw material and contemporary technology under the guidance of skilled professionals in compliance with the set quality standards. Also, our products are acknowledged for their features like robust construction, designs, longer life, rust proof nature, optimum quality, reliability, and durability. As per the requirements of the customers, these offered beauty and physiotherapy equipment can be customized.

Our staff follows a very systematic method enabling tender functioning of all the operations. Our infrastructural facility is segmented in more than a few departments and is well organized with the cutting edge science and developed machinery. We are into providing broad gamut of world classification products to our consumers which might be fine through all means. We thoroughly realise the intention with which these merchandise are bought and we put complete efforts to satisfy these wishes with whole devotion.

Spread over a wide area, we have large infrastructure facility that is known for carrying out standardized operations. In addition, the complete range is checked on various quality constraints before final deliverance at clients' premise. Known for our high production capacity, our infrastructure facility is upgraded on the regular basis in order to meet the various requirements of ever-changing market trends in the most efficient manner.

Primary Competitive Advantages

We are a customer centric organization and all our products and services are governed by stringent customer regulations and specifications. We work on various aspects of business and try to maintain a standard quality in overall business operation. Some of the salient features of our business operation are as follows:

- Experienced team of professionals
- Well-furnished infrastructural facility
- Total quality management
- Complete client satisfaction
- Ethical business policies
- Competitive pricing
- Immediate delivery

CONTENT

Co2 Fractional Laser Machine	4
Picosure Tattoo Removal Nd Yag Laser Machine	5
Q Switched Laser	
Desktop ND Yag Laser for tattoo removal + Black Doll	6
Desktop ND Yag Laser for tattoo removal + Black Doll	
1064nm+808nm+755nm Diode Laser Machine	7
810NM DIODE LASER HAIR REMOVAL SKIN REJUVENATION MACHINE ND YAG LASER	
808 Diode Laser Hair Removal	8
810 nm Diode Laser hair Removal Skin Rejuvenation Machine	
OPT Hair Removal Machine	9
E-LIGHT+SHR+OPT	
Desktop IPL Hair Removal Machine	10
Desktop Double OPT Handle	
Hair Removal Machine Desktop E-light Beauty Machine	
Lipo Hifu Slimming Machine 2017 Model	11
Fat Freeze Machine	
Double Handle Fat Freezing Machine	
Cryolipolysis OT-230	12
OPT/SHR & Eight hair Removal Machine Parameter	
Professional Cryolipolysis Slimming Machine	13
Portable Cryolipolysis Slimming Machine	
Body & Slimming Sculpture Beauty Machine	14
3D Lipo Slimming Machine	
5 In 1 Slimming Machine (Triangular Case)	15
Multipolar RF Fat Burst Slimming Machine	
Upgrade Version Medical Hifu	16
Upgrade Version Medical Hifu & Vaginal Tightening	
The Second Generation Luxury LED Beauty Machine	17
Mesotherapy Machines	18
Oxygen & BIO Skin Lifting Beauty Equipment	19
Oxygen & BIO Skin lifting Beauty Equipment	20
3 IN 1 Water Oxygen Skin rejuvenation Device	
Desktop 3 in 1 Water Dermabrasion Device	21
Mini Small Bubble Skin Rejuvenation machine	
Mini Small Bubble Skin Rejuvenation Machine	
Korean Small Bubble Moisture Machine	



Co2 Fractional Laser Machine

SI -201

Principle

CO2 fractional laser effect principle is 'focal photo-thermal effect'. Skin tissue has ability to repair itself. The skin was injured by trauma, but it could heal itself. The wound area is too large, the results of repair of scar formation. The wound is small, the skin can completely recover. Fractional laser every effect to skin only on one of the small region to generate thermal stimulation, retaining portion of skin is not affected, the role of tissue to the skin to stimulate collagen receipt hyperplasia and no effect of skin tissue can be used as thermal diffusion region, avoid the possibility of thermal injury and other side effects and can quickly promote skin healing, reach, and improve scars, skin firming pores purpose



Picosure Tattoo Removal Nd Yag Laser Machine

SI-202

Principle:

Pico Sure laser is using 775 wavelength laser with the collocation of honeycomb instantaneous lens (Focus), focus and collect the 70% laser energy, each laser provide 10% efficient shockwave function to the skin tissue, not only shatter the melanin cells of the epidermal layer, but also will conduct the shockwave to the deep cortex, start-up the repair of growth factor, promote renewal and proliferation of collagen and elastic fiber, accelerate skin melanin metabolism, accelerate the excretion of them, so as to solve the problem of deep pigment from the root.



Q Switched Laser

SI-203

Principle:

KEDUN-L6 laser skin treatment instrument treatment principle is the use of laser selective and Q-switched laser blasting principle, the precise does of specific wavelength of energy in a very short time alignment of a specific structure-color base: Dermal and epidermal ink, carbon particles and other exogenous pigment particles and endogenous pigment disease melanoma cells.

The pigment suddenly exploded by heat, was impacted for more small fragments, broken pieces of macrophage phagocytosis with macrophages into the lymph circulation and thus the pigment was excreted.



Desktop ND Yag Laser for tattoo removal + Black Doll

SI-204

Principle:

Laser emitted by the system has a strong penetration ability which allows it to reach the deep layer of dermis. The pigment particles absorb the tiny pieces, thus diminish the color density and get rid of it. So the appliance can effectively do away with mutant pigmentations and vascular tissue based on undamaged ambient tissue. This is called 'selective heat absorption' principle in medical field.



Desktop ND Yag Laser for tattoo removal + Black Doll

SI-205

Function:

1. OPT SHR Permanent hair removal;
2. Skin rejuvenation, wrinkles removal, and pigment removal;
3. Improving skin's situation: more smooth, tender, delicate and elastic;
4. Changing complexion dark: skin whitening and uniform color;
5. Removing the pigment and scar: removing the sun-induced freckle, speckle, etc.
6. Removing the acne and blain to imprint
7. Removing the vascular lesions and blood vessel



DIODE LASER HAIR REMOVAL SKIN REJUVENATION MACHINE

SI-206

wavelength	808nm	OE Mcustomized service
Working frequency	1~10hz	Maximum frequency 20Hz
Energy density	1-120j/cm2	1- 150j/cm2 optional
Handle piece cooling	TEC+sapphire	Sapphire size customizable
TEC cooling	12V 6A	TECpower customizable
Spot size	15cm x 15cm or 25 x 31mm	Spot size and window of light customizable
Output power	600w, 900w. 1200w	Laser power available
Laser bars	Germany laser Jenoptik, 150W/piece	
Pulse width	5ms~400ms	0~1000ms customizable
Output mode	Pulse output	Pulse output
Language	Standard English	customizable
display	8' color touch screen	10.4' optional
Protective temperature	55°C interface select	customizable
Water pump	DP-60 diaphragm rods 12V	imported pumpoptional
Flow sensor	At least 3.5L/min	
Water level switch	Test water level	
Cooling system	240 x 120cooling line(2pc) , 18 pcs channel	Compressorplan , TECplan
Cooling fan	4pcs fans , speed 3150r/min	12v fans
Case material	ABSshell/metal shell	customizable metal shell
Inter structure	White zinc hob/color zinc hob	customizable
Water tank volume	vertical3.6L/table type 2.5L	customizable
Environmental moisture	< 80%	< 80%
Environment temperature	Air conditioning room ' temperature ≤ 35°C	
voltage	220V	110Vavailable



Picosure Tattoo Removal Nd Yag Laser Machine

SI-207

Basic parameters:

755PicoSure Cellular Picosecond Laser

Laser Type: Picosecond Laser

Laser wavelength: 755,1064,532

Laser output:

Power supply: 220V / 110V

The highest energy: 3500mj / 2000mj / 1000mj

Operating frequency: 1-20HZ

Pulse width: 10ms

Power: 2000W

Cooling system: air-cooled

Spot diameter: 2-10mm

Net amount: 140KG



1064nm+808nm+755nm Diode Laser Machine

SI-208

Spot size: 13×13mm

Diode laser bar: 100W with 10 laser bars

Shots: 20 million shorts

Energy: 1-138J



810NM DIODE LASER HAIR REMOVAL SKIN REJUVENATION MACHINE ND YAG LASER

SI-209

Technical Parameter

Name	Multifunction yag laser machine
Display screen	8.4inch colour touch
Language on screen	seven different languages
Function	tattoo removal, remove lipline, eyeline, birthmark
and nevus spilus	
Laser wavelength	1064nm/532nm/1320nm(black moppet head)
Lamp longevity	800000shots
Power	800W
Energy	50-2000mJ (50 step)
Frequency	1-10Hz (1step)
Pulse width	1-100s
Control working or stop	foot switch
Cooling system	air+water circulation+semiconductor
Working environment	temperature: 18~23 /humidity: 10~90%
Voltage	110V/220V
Package size	130*65*61cm
Net weight	30kg
G.W	45kg
Package	wooden box with foam



810nm Diode Laser Hair Removal
Skin Rejuvenation Machine



Nd Yag Laser

808 Diode Laser Hair Removal

SI-210




810 nm Diode Laser hair Removal Skin Rejuvenation Machine

SI-211



810 handle



permanet hair removal



high power



in-motion design



Application

Remove lip hair, Beard shaved hair Back hair, arms hair,, Leg hair Bikini hair, Armpit, Chest hair Remove beard, Skin rejuvenation.

Advantage

1. High power 200W
2. At least 4000000 laser shots
3. energy output can reach 100 J/cm²
4. Fast 1-3 treatments and no pain
5. Permanent hair removal for all skin types and hair colors

Technical Specification

Screen: 10.4 inchinch
Laser type: diode laser
Number of bar: 10 bars
Frequency output: 1-5Hz
System: fast mode/normal mode for hair removal. Skin rejuvenation

Energy: 1-50J/cm²
Pulse: 10-200ms
Spot size: 12mm*12mm
Voltage: AC 110-220V 50Hz/60Hz
Cooling System: semi-conducotr + water circulation

OPT Hair Removal Machine

SI-212

Principle:

OPT perfect Sustainable constant light pulse technology to avoid energy attenuation, so that lesions quickly and accurately to achieve therapeutic energy, ensuring comfort and effectiveness of treatment, non-invasive hair removal, freckle, skin rejuvenation. once OPT treatment efficacy equivalent to 3-5 times of traditional therapy, painless and effective safer, shorter course of treatment.




E-LIGHT+SHR+OPT

SI-213

Principle:

1. Hair removal: painless fastly permanent hair removal painless, hairs through whole body (hairs in growing period, including tiny hairs with light color)
2. Freckle removing: removing freckles, chloasma, sunburn, age spots, acne marks and facial blemishes
3. Skin Rejuvenation: improving large pores, rough skin, tiny wrinkles, and restoring skin elasticity.
4. Telangiectasia treatment: redness, facial flush.
5. Improving dull complexion to whiten and uniform skin



Desktop IPL Hair Removal Machine

SI-214



1. Skin Rejuvenation: improving large pores, rough skin, tiny wrinkles, and restoring skin elasticity.
2. Freckle removing: removing freckles, chloasma, sunburn, age spots, acne marks and facial blemishes.
3. Hair removal: armpit hair, beard, hair limbs, bikini bob line etc effect the beauty hair

Desktop Double OPT Handle Hair Removal Machine

SI-215



1. Permanent Hair Reduction on all pigmented hair and all skin types-including tinny hair with light colour
2. Pain-free, Hair-Free - Comfortable Hair Removal.
3. Fast removing any kind and color hair permanently.
4. Skin rejuvenation and skin tightening.

Desktop E-light Beauty Machine

SI-216



1. Rejuvenation: improve the forehead, face, arms, abdomen, legs and other rough skin, enhance skin elasticity, whitening tender
2. Hair removal: permanent hair removal, hairs through whole body (hairs in growing period, including tiny hairs with light color)
3. Dispel spot: freckles, chloasma, freckle, chloasma, sunburn, senile plaque dermal patches, birthmarks and pigmentation
4. Dispel capillary vessel: congenital capillary vessel, capillary vessel after change skin, red face, erythema, skin allergies and sensitivities.
5. Acne removing: facial, shoulder, chest, waist, buttocks and other parts of bacterial acne



OPT/SHR & Eight Hair Removal Machine Parameter

SI-217

Item	Parameter
input voltage	110~220VAC±15%
current	10A
output power	1700W
Xenon lamp	50~100mm
Power size	depth 330xlenght 220xhighly 133mm
filters	590nm, 530nm, 640nm
operation screen	8.0ince colour touch
cooling system	semi-condctor, water circulation
spot size	15*50mm
cooling system	5 level
lamp longevity	One hundred and twenty thousand
package size	56*50&65
net weight	22kg



OPT/SHR & Eight Hair With double handle Removal Machine Parameter

SI-218

Treatment parameters:

Nursing care	The number of pulses	Pulse cycle	First pulse width	Secondary pulse width	interval	Light energy
Hair removal	3P-5P	1S-3S	Depth color: 1ms-3ms Light color: 3ms-5ms	Depth color: 1ms-3ms Light color: 3ms-5ms	20ms 25ms	Energy from 1 start up plus
Freckle removal	1P-2P	1S-3S	Depth color: 1ms-3ms Light color: 3ms-5ms	Depth color: 1ms-3ms Light color: 3ms-5ms	20ms 30ms	Energy from 1 start up plus
Skin rejuvenation	1P-3P	1S-3S	Partial white : 3ms-5ms Moderate : 2ms-5ms Partial black : 1ms-2ms	Partial white : 3ms-5ms Moderate: 2ms-5ms Partial black : 1ms-2ms	10ms 20ms 30ms	Energy from 1 start up plus
Wrinkle removal	1P-3P	1S-3S	Partial white : 3ms-5ms Moderate : 2ms-5ms Partial black : 1ms-2ms	Partial white : 3ms-5ms Moderate: 2ms-5ms Partial black : 1ms-2ms	10ms 20ms 30ms	
Acne removal	1P-3P	1S-3S	3ms-5ms 2ms-5ms 1ms-2ms	Partial white : 3ms-5ms Moderate: 2ms-5ms Partial black : 1ms-2ms	10ms 20ms 30ms	



OPT+Yag Laser Hair Removal Machine

SI-219

Technical Parameter

Wavelength	530-1200nm or 590-1200nm or 640-1200nm
Counter Handle	Yes.
Spot size	15*50j/cm2
Energy	1-50J/cm2
Continually shooting timer	1-30s
Pulser repetition rate	1-10 Hz
Pulse	Single
Pulse on	2-15ms
Water temperature tank	38 degree
Display system	8.4" colorful touch screen
Drive Power	2,000W
Power supply	220V, 50Hz or 110V, 60Hz
Laser wavelength	532nm/1064nm/Professional black moppet head
Laser pulse energy	1-2000mj
Laser plug frequency	1-10



360 Magneto Optic Hair Removal Machine

SI-220



Feature:

Display	: 8" colorful touch screen
Handle	: 1(HR)
Size	: 15*50nm2
Energy	: 1-50j/cm2
Interval time	: 1-4s
Pulse repetition rate	: 1-6 P/s
Plus time	: 1-9ms
Fast model plus	: 1-9ms
Frequency	: 1-6Hz
Energy	: 1-2000mj
Power	: 3000W
Packing size	: 70*58*125CM
Power supply	: 220V, 50Hz or 110V, 60Hz

360 Magneto Optic Hair Removal Machine

SI-221



Display	: 8" colourful touch screen
Wavelength	: 530/590/640nm
Size	: 15*50nm2/15*32nm2(random)
Energy	: 1-50j/cm2
Elight discharge time	: 1-4s
Elight Pulse repetition rate	: 1-6 P/s
Elight plus	: 1-9ms
OPT SHR plus	: 1-8ms
OPT SHR frequency	: 1-5s
Power	: 2000W
Packing size	: 119*61*59cm
Continuous working time	: 8house
Power supply	: 220V, 50Hz

Cryolipolysis OT-230

SI-222

Specification

Cryolipolysis temperature	-10 to 5 degree (controllable)
Heat temperature	37°-45°
Vacuum Power	0-100KPa
Radio Frequency	5Mhz high frequency
LED wavelength	650nm
Cavitation frequency	40Khz
Cavitation modes	4 types of pulse
Lipo Laser length	650nm
Energy controlled	Yes
Laser modes	AUTO,M1,M2,M3
Machine Display	8.4 inch touch screen
Handle Display	3.5 inch touch screen
Cooling System	Semiconductor +water+air
Input voltage	220~240V/100-120V, 60Hz/50Hz
Packing size	123*60*59CM



Double Cryolipolysis Machine

SI-223

Speci·cation

Voltage	AC110V/220V 50-60Hz
Power	< 400w
Temperature	-14 °C~14° C
Pressure output	0-900Kpa
Handle working size	5*13cm or 7*16cm
Screen	8.4inch touch screen
Body RF	1Mhz
Cryolipolysis	-14 °C~5° C
Cavitation frequency	40Khz
Net weight	20KG
Gross weight	28KG



Lipo Hifu Slimming Machine 2017 Model

SI-224

Theory:

High Intensity Focused Ultrasound penetrates focused and volumetric ultrasonic energy at precise depth (13mm) in subcutaneous fat layers targeting fat tissue causing a non-thermal effect. The surrounding structures such as the epidermis, blood vessels and nerves will not be damaged. Meanwhile the metabolized contents (triglyceride, fatty acid) of cells will excrete outside body auto-matically through blood circulation and lymphatic drainage without harm to the human body.

Treatment scope:

Treatment area has an inch or 2.5cm of fat to treat, abdomen, waist, buttock and both sides of the thighs site.



Handpiece



0.8cm Emitter



1.3cm Emitter



Fat Freeze Machine

SI-225

Theory:

1. Minus 20' C--10' c freezing vacuum sculptor
2. 40KHZ strong wave cavitation
3. RF fat dissolve



Double Handle Fat Freezing Machine

SI-226

Theory:

A collection of "frozen melt fat", negation+light carve" "fat cavitation" three powerful technologies in one, in-depth skin fat cells under the intracellular effects of vibration and produce a hole, which would makethe fat cells resonance burst into positioning, dissolving fat and reducing the fat of orange pel



Professional Cryolipolysis Slimming Machine

SI-227

Specification

Model name	4S
Cryolipolysis temperature	5° to -10° adjustable
Heat temperature	About 40centi-degree
Vacuum	-100KPa
Vacuum mode cellulite	4 modes for different
PDT wavelength	650nm
Machine Display	10.4 inch touch screen
Handle Display	3.5 inch touch screen
Power supply	110V-120V/220V-240V, 50/60Hz
Gross Weight	70KG
Packing size (wooden box)	118*66*67cm



Portable Cryolipolysis Slimming Machine

SI-228

Specification

Cryolipolysis temperature	: 5° to -10° adjustable
Vacuum	: 102KPa
PDT wavelength	: 650nm
Display	: 8.4 inch touch screen
Powersupply	: 110V-120V/220V-240V, 50/60Hz



Body & Slimming Sculpture Beauty Machine

SI-229

Wave fat system

With collective strong sound wave head, strong sound wave of 40KHZ may be emitted to vibrate fat cells at top speed and produce numerous vacuum air pockets inside and outside the fat cells, robustly impact fat cells to generate introverted blast and disintegrate triglyceride into glycerol and free fatty acids. Then RF waves at frequency of 5MHZ are used for exhausting the integrate glycerol and free fatty acids through hepatoenteral circulation. Finally, vacuum RF and energy electrode are used for positioning and tightening fat. It physics, it is known as "cavitation". Micropore introverted blast inside and outside cell my lead to enhanced molecular motion and a higher energy level and this will finally cause fat cell rupture and thereby achieve the effects of body building and losing weight.



3D Lipo Slimming Machine

SI-230



Technical parameter

40KHZ Ultrasonic Head Frequency	40KHZ
Vacuum RF Head Frequency	1MHZ
Quadrupole RF Head Frequency	1MZH
BIO Head RF Frequency	1MZH
bi~polar face RF Head Frequency	1MHZ
Power	350W
Voltage	110AVC~240VAC
Working Temperature	10°C~35°C
G.W.	35KG
Meas	117*59*64cm





5 In 1 Slimming Machine (Triangular Case)

SI-231

Function

- 1 Face and Neck Lifting & Remodeling
- 2 Tightening and lifting flaccid areas of the face
- 3 Reduces wrinkles/left sagging Skin
- 4 Enhances luminosity of the skin
- 5 Fat Reduction
- 6 Body remodeling
- 7 Cellulite Treatment and skin partial reconstruction
- 8 Increase CollagenLevel
- 9 Tightening looseskin



Multipolar RF Fat Burst Slimming Machine

SI-232

Function

1. The facial wrinkle: the forehead wrinkles, glabella wrinkles, crow's feet, nose fissured lips, ears ago wrinkles, lower eyelid wrinkles and neck wrinkles.
2. promote and tight the eye, forehead, cheeks, neck lift solid; shrink pores, remove the double chine
3. body shaping: slimming shaping, eliminate excess fat, skin, fade striae of pregnancy.
4. Laser fat soluble fat: laser fat soluble fat, smooth skin without leaving marks, safety and no pain



Upgrade Version Medical Hifu

SI-233

Theory:

High Intensity Focused Ultrasound (HIFU) directly delivers heat energy to skin and subcutaneous tissue that can stimulate and renew the skin's collagen and thus consequently improving the texture and reducing sagging of the skin. It literally achieves the results of a facelift or a body lifts without any invasive surgery or injections, moreover, an added bonus of this procedure is that there is no downtime. and also, it works equally well for people of all skin colours, in contrast to that of lasers and the intense pulse lights




- 1.5mm
- 1.5mm
- 3.0mm
- 4.5mm
- 13mm



Upgrade Version Medical Hifu & Vaginal Tightening

SI-234



HOT HIFU MACHINE (MEDICAL VERSION)



The Second Generation Luxury LED Leauty Machine

SI-235

Theory:

Improve slack and rough skin; whiten and rejuvenate skin, moisten face. Despel freckle, sunspot, pigmenga-tion, lower skin damage from daylight. Shrink rough pore,. Lighten red mark and capillary vessel enlarge-ment of the erythema acne. Smoothen tiny wrin-kle, tighten slack and aging skin. Treat whelk, mark of acne, hair follicle inflammation. Skin doom, no shine for poor metabolism, circulation, Rrcover the face nerve with light pathological changes.



2017 X7 Led & Skin Care Machine



Led+skin Analyzer + Face Care Machine



New Multifunction Hair Regrowth Machine



Mesotherapy Machines

SI-236



Needle Free Mesotherapy Face Machine



Technical Specification
Diamond Dermabrasion Probe



Meso Therapy Machine



Facelift Probe



Negative Bar



Cooling Probe



Vacuum Detox Probe



Electroporation Probe



Under eye care probe



Radio Frequency Probe



Physiotherapy Shockwave Therapy



Portable Cryolipolysis



Latest Bubble Beauty Machine



Mini Omega Injection Sprayer Machine



Body Analysis System

Oxygen & BIO Skin Lifting Beauty Equipment

SI-237



Function:

SPA water dermabrasion machine completely changed the traditional way, which is hand-cleaning skin relying on individual's practice skills, SPA Skin beauty machine, using vacuum suction mode controlled by intelligent process, through the combination of products and equipment, deep cleaning the skin and pores off horny, acne, blackheads and other impurities in very short period of time, And improve the deep absorption of nutrition products, promote pores tightening, smooth skin, increase skin moisture, and make your skin whitening, moisturizing and good texture.

Application:

1. Improve symptoms such as aging and loosening skin, big pore, slender wrinkles.
2. Improve pigmentary pathological changes, such as freckle, sunburn, senile plaques.
3. Improve dark complexion caused bad metabolism or poor circulation.
4. Repair and nurse the damaged skin.
5. Do detumescence, inflammation diminishing and scare elimination treatment for all acne effectively




Oxygen & BIO Skin lifting Beauty Equipment

SI-238



Function:

1. Deep cleaning, Oily skin improvent
2. Scare removal: All kinds of scars suchas scar leftby laser, burn and surgery etc.
3. Acne: blain acne, scabby acne allergic acne, papilla acne, improving the appearance of lipidic skin and acne pit.
4. Skin care: skin whitening and softening, facial lifting and tightening, eye bad and black eye circle removing, weary skin and gloom yellow skin improving
5. Wrinkle reduction: reduce wrinkles around canthus furrow.
6. Hair regrowth: has good effect for Alopecia areata, baldness and loss of hair etc
7. To improve allergic skin
8. Water replenishing to skin.



3 IN 1 Water Oxygen Skin rejuvenation Device

SI-239

Function:

1. Spot elimination: Good for removing scald, scare left after surgery, acne makes and pockmarks;
2. Acne elimination: Good for treatment on acne, inflammation, pimple, and blackead.
3. Clean and provide water: Clean from the deep part of the skin so that make the skin cutin layer, minimally invasive invasive scare, blackhed remove, at the same time supply theskin with sufficient water;
4. Tenderize skin: Helpful to wipe to wipe off the patch, pigmentation and good for pore contraction so as to whiten skin and make the skin more elastic and shining. The condition of the skin all over the body would be improve.
5. Smooth wrinkles: Beneficial to smooth wrinkles around lips and eyes, wrinkles between eyebrows and smooth striae gravidarum.
6. Hair growth: Good for resist seborrheic alopecia, alopecia areata and have a effective caring onhair;
7. Bring a young state of the skin: Suitable to deal with pimple, facial holes, fat granules, acne marks, and various patches.



Desktop 3 in 1 Water Dermabrasion Device

SI-240



Function

1. Ideal for normal skin
2. Ideal for sensitive skin, acne pimple Etc.
3. Injecting nutrient liquor products
4. Facial cleansing
5. Cleaning with nutrient liquor products



Mini Small Bubble Skin Rejuvenation Machine

SI-242



Function

1. acne, folliculities, seborrheic alopecia, clear skin mites and allergicsource;
2. skin whitening improve dark and yellow skin; improve skin texture;
3. deeply clean the skin, moisturize skin nourishing, oxygenating;
4. to wipe out balckhead; improve the skin slack, minimise the pores, increase the transparency of the skin;
5. preoperative and postoperative nursing care of patients with skin reconstruction and non stripping skin reconstruction

Hydra Facial Machine

SI-241



Function

1. acne, seborrheic alopecia, folliculites, mites clear, clear skin allergens;
2. skin whitening, improve skin dull, yellowish, improve skin texture;
3. Deep clean the skin, moisturizing skin, nourishing;
4. Julep, improve loose skin, tighten pores, increase skin transparency;
5. The preoperation and postoperation care for ablative skin reconstruction and non-ablative skin reconstruction;
6. shaping firming skin, shrink pores, improve double chin;

Korean Small Bubble Moisture Machine

SI-243



Function

1. acne, seborrheic alopecia, folliculites, mites clear, clear skin allergens;
2. skin whitening, improve skin dull, yellowish, improve skin texture;
3. Deep clean the skin, moisturizing skin, nourishing;
4. Julep, improve loose skin, tighten pores, increase skin transparency;
5. The preoperation and postoperation care for ablative skin reconstruction and non-ablative skin reconstruction;
6. shaping firming skin, shrink pores, improve double chin;



SURJEET INTERNATIONAL



Surjeetinternational81@gmail.com

9911442933,7529979797

H.O.



522 Block M, Guru Harkishan Nagar,
Paschim Vihar, New Delhi-110087

B.O.



www.pryogaphysionslim.com

